Urine Containers

With screw cap closure



Urine containers from PE

- ☐ Brown translucent material light protection
- Volume graduation
- ☐ Self-adhesive label enclosed (2.0 I container) detailing instructions for use and sample identification (see illustration)
- ☐ Large aperture
- ☐ Practical moulded carrying rim
- ☐ Clean room manufacture, cap assembly immediately after production, thus no falsification of analytical values by foreign particles
- ☐ Container is coloured with heavy metal-free pigment

Urine containers from PE with assembled yellow cap

Order Number	Volume	Description	Packaging Pcs./Case
77.593	2.0 Litre	self adhesive label enclosed	30
77.576	2.5 Litre	without label	30

ml hours illet in the morning:
hours
illet in the morning:
this collection container upty bladder into the it the same hour as stainer from scale and ided space.

Self adhesive label

detailing instructions for use and sample identification (see illustration)

Label dimension: 115 x 75 mm

Order Number: 84.1398



Urine containers from PE with screw cap

Order Number	Volume	Description	Packaging Pcs./Case
77.582	0.5 Litre	brown translucent material, light protection	100
77.577	0.25 Litre	white	120